



CONSUMABLES | CLEANERS | NEUTRALIZERS

Formalin Neutralizer Powder

Mopec Formalin Neutralizer Powder is an EPA validated system for the disposal of solutions containing formaldehyde and glutaraldehyde (up to 4% w/v). Mopec Neutralizer utilizes fast acting reagents that neutralize aldehydes quickly and effectively. Mopec's proprietary aldehyde reduction system achieves required pH levels through a novel approach of buffers and stabilizers, making an approved neutralized solution the first time, every time.

SPECIFICATIONS

Instructions for Use

1. Collect waste in a designated collection container.
2. When collecting waste, be sure to use a funnel with filtering device to remove tissue debris.
3. Add the entire contents of one NEUTRALIZER packet to each gallon (128-fl.oz) or to 4-liters of aldehyde waste.
4. Secure lid on designated collection container and agitate container to thoroughly mix.
ENSURE THE ENTIRE AMOUNT OF NEUTRALIZER HAS FULLY DISSOLVED.
5. Allow the mixed solution to stand for 15-20 minutes to complete the neutralization process.
6. Test treated solution for any residual aldehyde content and proper pH levels prior to disposal.
7. Upon passing aldehyde content and pH tests, the treated waste is now ready to be disposed down the drain. Consult with local and regional authorities to ensure complete and accurate classification of hazardous waste.
8. After disposal, rinse collection container with cold tap water to prepare for next use.

Results

The treated waste is safe to dispose if the following:

- pH range: 6-9
- Aldehyde Test: Negative

Packaging

| SKU | Volume | Packet Dimension | Case Dimensions |
|----------|-----------------|----------------------|------------------------|
| CDN15FNP | 30 Packets/Case | 9"L x 5.75"W x 1.5"H | 12.5"L x 10.5"W x 11"H |



www.mopec.com
800.362.8491
21750 Coolidge Hwy.
Oak Park, MI 48237