



Surgical Lights

BB108-BB113

* Specifications of modified models may vary from the model shown.

Specifications

Bulbs & Brightness

- 86,000 lux (7990 fc) at 24" (61 cm)
- 3300 K color temperature
- CRI (Color Rendering Index) of 94
- Light uses three 50-watt halogen bulbs, each with 2000-hour average bulb life.
- Bulbs are wired in parallel to ensure continuous operation

Mechanisms

- Focuses with central removable, autoclavable SteriHandle
- Friction knob allows precise tool-free adjustment and positioning of the articulating arm
- Cool operation ensured with heat filters and dichroic coatings
- The light head's unique shape helps reduce shadows on examination field

- Design allows removal or replacement of SteriHandle with one hand (for ensuring sterile procedures)
- Floorstand units include footswitch, locking casters and cord wrap

Other Standard Features

- Designed to meet California Seismic Codes
- 120V model standard; 230 / 240V models also available
- UL 60601-1, EN 60601-1-2, IEC 60601-2-41 and CAN/CSA C22.2 NO. 601.1 M90
- Assembled in USA



MODEL #	DESCRIPTION
BB108	Floorstand
BB109	Double Ceiling Mount
BB110	Single Ceiling Mount
BB111	Wall Mount
BB112	Single Fastrac Mount
BB113	Double Fastrac Mount