



# Autopsy Cart Chassis

## DA000

\* Specifications of modified models may vary from the model shown.

### Specifications

#### Standard Features

- Chassis only (cart top not included)
- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Four heavy duty 8" casters, 3 total lock and 1 directional lock

#### Performance

- Weight capacity: Approximately 500 lbs

#### Dimensions

- 53" L x 31.25" W  
36.7" H - 34.7" H cart top post height to provide slope

#### Compatibility

- Autopsy Sinks: CA, CB and CC Series
- Cart Tops: GA100 and GA300

