



Autopsy Cart

DA100

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- Removable stainless steel cart top (GA100)
- Includes drain hole and plug
- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Four heavy duty 8" casters, 3 total lock and 1 directional lock

Performance

- Weight capacity:
Approximately 500 lbs

Dimensions

- 80.5" L x 32" W
37.4" H - 34.3" H Sloped top height

Compatibility

- Autopsy Sinks: CA, CB and CC Series
- False Frame & Cover: GA110

