



## Autopsy Cart with Perimeter Trough

DA300

\* Specifications of modified models may vary from the model shown.

### Specifications

#### Standard Features

- Removable stainless steel cart top (GA300)
- Perimeter trough on cart top
- Includes drain hole and plug
- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Four heavy duty 8" casters, 3 total lock and 1 directional lock

#### Performance

- Weight capacity:  
Approximately 500 lbs

#### Dimensions

- 80.5" L x 32" W  
37.4" H - 34.3" H Sloped top height

#### Compatibility

- Autopsy Sinks: CA, CB and CC Series
- False Frame & Cover: GA110

