



Tilting Autopsy Cart Chassis

DB000

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- Chassis only (cart top not included)
- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Four heavy duty 5" casters, 3 total lock and 1 directional lock

Performance

- Weight capacity:
Approximately 500 lbs

Dimensions

- 53" L x 31.25" W
fixed end height of 35"

Range of Motion

- Tilting cart end lifts from
36" - 42"

Controls

- Hand crank with hydraulic lift system for tilting mechanism

Compatibility

- Autopsy Sinks: CA, CB and CC Series
- Cart Tops: GA100 and GA300

