



Tilting Autopsy Cart with with Perimeter Trough

DB300

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Removable cart top (GA300)
- Perimeter trough on cart top
- Includes drain hole and plug
- Four heavy duty 5" casters, 3 total lock and 1 directional lock

Performance

- Weight capacity: Approximately 500 lbs

Dimensions

- 80.5" L x 32" W
fixed end height of 34.88"

Range of Motion

- Tilting cart end lifts from 36.44" - 42.44"

Controls

- Hand crank with hydraulic lift system for tilting mechanism

Compatibility

- Autopsy Sinks: CA, CB and CC Series
- False Frame and Cover: GA110

