



Bariatric Autopsy Cart

DC400

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- Non-removable bariatric top
- All welded stainless steel frame, type 304 1.5" SQ. 11 gauge
- Includes drain hole and plug
- Four heavy duty 8" casters with 3 total lock and 1 directional lock

Performance

- Weight capacity:
Approximately 1000 lbs

Dimensions

- 80.5" L x 40" W
- Sloped top height:
37.4" - 34.3"

Compatibility

- Autopsy Sinks: CA, CB and CC Series
- False Frame & Cover: DC410

