



# Tilting Bariatric Autopsy Cart

DC500

\* Specifications of modified models may vary from the model shown.

## Specifications

### Standard Features

- Non-removable bariatric top
- All welded stainless steel frame, type 304 1.5" SQ. 11 gauge
- Includes drain hole and plug
- Four heavy duty 8" casters with 3 total lock and 1 directional lock

### Performance

- Weight capacity:  
Approximately 1000 lbs

### Dimensions

- 80.5" L x 40" W
- Fixed end height: 34"

### Range of Motion

- Tilting cart end lifts from  
35.7" - 41.7"

### Compatibility

- Autopsy Sinks: CA, CB and CC Series
- False Frame & Cover: DC410

