



Saddle Chassis Autopsy Cart for 27" Trays

DD050

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Four heavy duty 5" casters with 3 total lock and 1 directional lock

Performance

- Weight capacity:
Approximately 375 lbs

Dimensions

- 79" L x 28.75" W
- Sloped top height:
35" - 34"

Compatibility

- Autopsy Sinks: CA, CB and CC Series
- 27" Body Tray: JC101
- Forklifts: JD710, JD305, JD405

