



## Roller Autopsy Chassis for 23" Tray

DE000

\* Specifications of modified models may vary from the model shown.

### Specifications

#### Standard Features

- All welded stainless steel frame, type 304 1.5" SQ. 16 gauge
- Chassis only; cart top/tray not included
- Conveyor roller transfer style top
- Four heavy duty 8" casters with 3 total lock and 1 directional lock

#### Performance

- Weight capacity:  
Approximately 375 lbs

#### Dimensions

- 79" L x 27.25" W
- Sloped top height:  
35" - 32.6"

#### Compatibility

- Autopsy Sinks: CA, CB and CC Series
- 23" Body Tray: JC100

