



Creased Top Dissection Table

HA200

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel construction
- Creased non-removable top
- Built-in tapered slope
- Includes drain hole PVC drain valve
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

Performance

- Weight capacity:
Approximately 350 lbs

Dimensions

- 78" L x 27" W
- Sloped top height: 36" - 35"

Optional Features

- Post style leg rests (HO005)
- Post style book holder (HO006)
- Stainless steel bucket (BA005)

