



## Dissection Cart with Removable Tray

HA300

\* Specifications of modified models may vary from the model shown.

### Specifications

#### Standard Features

- All welded stainless steel frame; type 304 1.5" SQ. 16 gauge
- Conveyor roller transfer top
- Removable 27" body tray (JC101)
- Includes drain hole PVC drain valve
- Built-in slope
- Lower shelf
- Includes 2.5 gallon carboy and carrier
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

#### Performance

- Weight capacity:  
Approximately 350 lbs

#### Dimensions

- 80" L x 30.2" W
- Sloped top height:  
35.8" - 34.4"

#### Optional Features

- Post style book holder (HO006)

