



# Hinged Cover Dissection Table

HB100

\* Specifications of modified models may vary from the model shown.

## Specifications

### Standard Features

- All stainless steel construction
- Recessed non-removable top
- Hinged cover with swing latches to hold open when in use
- Built-in tapered slope
- Includes drain hole PVC drain valve
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

### Performance

- Weight capacity:  
Approximately 350 lbs

### Dimensions

- 86" L x 30" W
- Sloped top height: 36" - 35"

### Optional Features

- Bridge style leg rests (HO001)
- Bridge style book stand (HO002)
- Bucket holder (HO003)
- Post style leg rests (HO005)
- Post style book holder (HO006)
- Fixed stainless steel shelf (HO024)
- Stainless steel bucket (BA005)

