



Downdraft with Tilt Elevation Dissection Table **HB300**

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel construction (cart top only)
- (4) removable perforated grid plates and a recessed top
- Twin pedestal hydraulics for easy foot activated elevation
- 6" dia. vent stub with damper
- PVC drain valve
- Four heavy duty 8" casters with total lock

Performance

- Downdraft system removes formalin vapors down and away from user
- Weight capacity:
Approximately 500 lbs

Dimensions

- 86" L x 30" W
- Elevates from 31.25-44.75" H

Power

- Exhaust requirements:
350-500 CFM

Controls

- Brake lever on both sides and one on end
- Elevating pedals on both sides

Optional Features

- Bridge style leg rests (HO001)
- Bridge style book stand (HO002)

