

Elevating Downdraft with Hinged Top Dissection Table


HB450

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel construction
- Removable perforated grid plate and a recessed top
- Hinged cover with swing latches to hold open when in use
- 6" dia. vent stub with damper
- Includes drain hole PVC drain valve
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included



Performance

- Downdraft system removes formalin vapors down and away from user
- Weight capacity: Approximately 350 lbs

Dimensions

- 86" L x 30" W
- Sloped top height (min): 32-31" H
- Sloped top height (max): 43-42" H



Power

- Exhaust requirements: 220-340 CFM
- Static pressure: .34 - .67" WG

Controls

- Hand crank elevation control

Optional Features

- Bridge style leg rests (HO001)
- Bridge style book stand (HO002)
- Post style leg rests (HO005)
- Post style book holder (HO006)
- Stainless steel instrument tray holder (HO030)