



Two Tier Portable Cadaver Rack For 23" Trays

IF23201

* Specifications of modified models may vary from the model shown.

Specifications

Construction

- All stainless steel construction
- Polymer roller assembly
- Four heavy duty 5" casters, 3 total lock and 1 directional lock
- Rear bars with bumpers for loading the body tray
- Tray latch at each tier

Dimensions

- 81.5" L x 27.25" W x 44.2" H
- 14.5" body clearance between trays

Configuration

- 23" tray x 2 tiers high x 1 bay

Compatibility

- Accommodates 23" Body Trays:
JC023 (sold separately)
- Works with Cadaver Lifts:
JD301, JD501, JD5000 and JD5001

