



Hydraulic Stainless Steel Top Cadaver Carrier with False Frame & Blue Cover **JA614**

* specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel fixed top
- Twin pedestal hydraulics for easy foot activated elevation
- Four heavy duty 8" casters with total lock
- Bumper guards
- Removable aluminum false frame with heavy duty blue opaque cover (JZ007)

Performance

- Weight capacity:
Approximately 700 lbs

Dimensions

- 81" L x 30" W

Range of Motion

- Elevates from 21.25" to 34" (carrier only)
- Elevates from 37.25" to 50" (with false frame)

Controls

- Brake lever on both sides and one on end
- Elevating pedals on both sides

