



Better By Design™

Necropsy



Equipment and Accessories for Animal Research

ANATOMIC PATHOLOGY | MEDICAL EXAMINER/MORGUE | ANATOMY LABS | NECROPSY

Enhance Your Research With Innovative Technology

TRIMMING WORKSTATIONS

Our trimming tables and workstations feature fully welded construction with heavy-duty undercarriage subframes. Various style table tops, ventilation, rotation and elevation options are available.



TRANSPORTATION

We support specimen handling efforts with a full selection of racks, trays, carriers, pallet systems and lifts designed to work in sync with each other.



REFRIGERATION

Our cooler selection consists of telescopic, conveyor, walk-in, small door, and specimen locker with modularity options available. Select coolers include Sanisteel®, an antibacterial coating embedded into the walls of select models.



CASEWORK

Our casework selection includes wall mount, base, corner and specialty cabinetry and may be used to store specimens or accessories.



FUME HANDLING

Mopec has a multitude of fume handling solutions available. We offer a variety of hoods, adapters, modules, workstations, storage and work cells to suit your necropsy needs. Our potassium permanganate and charcoal filters each have the ability to absorb over 50 different chemicals.



ACCESSORIES

Hygiene and instrumentation are essential to the necropsy process. Mopec offers a variety of cleaners, disinfectants, deodorizers and specialty instruments to complement your lab's equipment.



DESIGN & ENGINEERING FOR THE ANIMAL RESEARCH ENVIRONMENT

Mopec engineers will work as a partner with your project team to create or renovate your necropsy lab. We provide ideas and concepts illustrating the facility workflow with full 3D rendered models. Our team will provide detailed product information to the end user, facility engineer and contractors. We also provide on-site project assistance and review.



PLANNING & WORKFLOW

- NEEDS ANALYSIS
- PROPOSAL REVIEW
- SYSTEMS APPROACH
- DESIGN
- VALIDATION
- DELIVERY
- INSTALLATION
- POST SALES



FABRICATION FACILITY

- Fully computerized fabrication equipment
- Live job tracking through the fabrication process
- Lean manufacturing is used to streamline workflow
- ISO Certified 9001



CUSTOM ENGINEERED AND MANUFACTURED TO MEET YOUR REQUIREMENTS

- Custom equipment built to your needs and standards
- Conceptual product design
- Fully engineered plans
- Value engineering
- Ergonomics and safety
- Products engineered to solve problems


Mopec's turnkey systems approach to necropsy lab development creates a comprehensive environment that focuses on safety, ergonomics, versatility and work flow. Our customized equipment is designed to accommodate large and small animal necropsies.

About Mopec



Better By Design™

Mopec is committed to designing and manufacturing equipment and accessories to enhance your facility's safety, flexibility and productivity. Our technical prowess stems from true American craftsmanship and high touch service. Mopec is recognized globally for giving clients an unrivaled experience including design, manufacturing, installation and post-sales support. The industries we serve include anatomic pathology, medical examiner/morgue, anatomy labs and necropsy.



INDUSTRIES WE SERVE



Anatomic Pathology



Medical Examiner/Morgue



Anatomy Labs



Necropsy

VISIT WWW.MOPEC.COM FOR MORE INFORMATION

21750 Coolidge Highway • Oak Park, MI 48237 USA
+1 800-362-8491 • 248-291-2040 • Email: info@mopec.com

Equipment shown may include optional features. Please visit our website for more information.

Copyright ©Mopec 2015