



Creased Top Dissection Table with Drawers

HA400

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel construction
- Non-removable top with a 1" recessed vertical edge
- Includes drain hole PVC drain valve
- Built-in tapered slope
- Three stainless steel drawers
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

Performance

- Weight capacity:
Approximately 350 lbs

Dimensions

- 87" L x 30." W
- Sloped top height:
36" - 34.5"

Optional Features

- Bridge style leg rests (HO001)
- Bridge style book stand (HO002)
- Bucket holder (HO003)
- Post style leg rests (HO005)
- Post style book holder (HO006)
- Fixed stainless steel shelf (HO024)
- Stainless steel bucket (BA005)

