



# Economy Dissection Table

## HA210

\* Specifications of modified models may vary from the model shown.

## Specifications

### Standard Features

- All stainless steel construction
- Creased non-removable top
- Includes drain hole PVC drain valve
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

### Performance

- Weight capacity:  
Approximately 350 lbs

### Dimensions

- 80" L x 30" W x 36" H

### Optional Features

- Bucket Holder (HO003)
- Post style leg rests (HO005)
- Post style book holder (HO006)
- Stainless steel bucket (BA005)

