



Economy Dissection Table

HA210

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- All stainless steel construction
- Creased non-removable top
- Includes drain hole PVC drain valve
- Four heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

Performance

- Weight capacity:
Approximately 350 lbs

Dimensions

- 80" L x 30" W x 36" H

Optional Features

- Bucket Holder (HO003)
- Post style leg rests (HO005)
- Post style book holder (HO006)
- Stainless steel bucket (BA005)

