

Chalk Calvaria Elevator



Designed with
function and
comfort in mind.

Features

The Chalk Calvaria Elevator's penetration edge is only 10mm long, which prevents unintentional damage of the cranial cavity and laceration to the brain.

The heavy-duty strike plate and handle are mounted at a right angle to the penetrating edge. Therefore making the use of a mallet easier and allowing the user to apply torque when removing the calvarium.

With the Chalk Calvaria Elevator any discomfort from the torque method is eliminated, which in turn will reduce user injury.

*Patent Pending

