



Grossing Station Imaging System With Stereoscopic Viewing

Specifications

Standard Features

- 20 MP Quality Imaging
- Auto Focus
- 10x Optical Zoom
- 12x Digital Zoom
- Digital Output: Frames/Sec
1920 x 1080p @ 25,30,50,60
1920 x 1080i @ 50,60
1280 x 720p @ 25,30,50,60
- Max Fof V 28" x 36" and
Min Fof V .185" x 1.85"

Dimensions

- **PathCam**: 4.50" length,
2.25" diameter & approx. 375 g
- Blue cutting board: 18.5" x 27"

Performance

- Mounted in hood and/or
tethered and retractable
- Measure distances and
area in Snapshots
- Grid tool: size biopsies and
map structure location
- Create system folders/files
- Record independent voice
narration
- Record video files with or
without voice narration

Power

- 12V DC
- (7V to 15VDC)
- -20C ~ 60C

Compatibility

- Windows 10 open architecture
- Barcode reader compatible
- Folders/files transferable to any
OS system or network
- Works with any Windows word
processor or spreadsheet App.
- Snapshots archived as JPEG's
- Voice narrations archived as WAV's
- Video recordings archived as AVI's
- Videoconference compatible with
3rd party webinar software

Optional Features

- Retractable with **PathCam** control on
handle: four (4) programmable buttons
- Foot Pedal **PathCam** control:
four (4) programmable pedals
- Voice Command **PathCam** control
- Washable keyboard and mouse
- 21" 4K Touch Screen Monitor
- 65" 4K **Ultra-D™** (Glasses-Free 3D)
- PC with Windows 10: 8 USB,
1 HDMI, and 1 VGA port
(required and available for purchase)

