



## Grossing Station Imaging System

### With Stereoscopic Viewing


### Specifications

#### Standard Features

- 20 MP Quality Imaging
- Auto Focus
- 10x Optical Zoom
- 12x Digital Zoom
- Digital Output: Frames/Sec  
1920 x 1080p @ 25,30,50,60  
1920 x 1080i @ 50,60  
1280 x 720p @ 25,30,50,60
- Max Fof V 28" x 36" and  
Min Fof V .185" x 1.85"

#### Dimensions

- *PathCam*: 4.50" length,  
2.25" diameter & approx. 375 g
- Blue cutting board: 18.5" x 27"



Independent voice recording

#### Performance

- Mounted in hood and/or  
tethered and retractable
- Measure distances and  
area in Snapshots
- Grid tool: size biopsies and  
map structure location
- Create system folders/files
- Record independent voice  
narration
- Record video files with or  
without voice narration

#### Power

- 12V DC
- (7V to 15VDC)
- -20C ~ 60C



*PathCam* field of range

#### Compatibility

- Windows 10 open architecture
- Barcode reader compatible
- Folders/files transferable to any  
OS system or network
- Works with any Windows word  
processor or spreadsheet App.
- Snapshots archived as JPEG's
- Voice narrations archived as WAV's
- Video recordings archived as AVI's
- Videoconference compatible with  
3rd party webinar software

#### Optional Features

- Detachable with *PathCam* control on  
handle: four (4) programmable buttons
- Foot Pedal *PathCam* control:  
four (4) programmable pedals
- Voice Command *PathCam* control
- Washable keyboard and mouse
- 21" 4K Touch Screen Monitor
- 65" 4K Ultra-D™ (Glasses-Free 3D)
- PC with Windows 10: 1 on-board USB3,  
7 USB2, 1 HDMI, and 1 VGA port  
(required and available for purchase)