



# Covered Dissection Table with Lid Assist

## HB Series

\* Specifications of modified models may vary from the model shown.

### Specifications

#### Standard Features

- Hinged covers with lid assist to safely and easily open and close
- All stainless steel construction
- Recessed, non-removable top
- Level cart surface with sloped recess
- Includes drain with PVC ball valve and drain tube
- Four heavy duty 8" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

#### Performance

- Hoods cover cadaver approximately up to 250 lbs.
- Weight capacity: Approximately 350 lbs

#### Optional Features

- Bridge style leg rests
- Bridge style book stand
- Removable bucket holder
- Stainless steel bucket
- Post style leg rests
- Post style book holder
- Removable recessed 48" shelf
- Lower fixed drawers



|         |           | Features             |           |                          |
|---------|-----------|----------------------|-----------|--------------------------|
|         | Model No. | Ventilation          | Elevation | Dimensions               |
| Product | HB1100    |                      |           | 86"L x 32.5"W x 36"H     |
|         | HB1410    | 6" Single Connection |           | 86"L x 32.5"W x 36"H     |
|         | HB1420    | 6" Dual Connection   |           | 86"L x 32.5"W x 36"H     |
|         | HB1510    | 6" Single Connection | 32"-38"   | 92"L x 32.5"W x 32"-38"H |
|         | HB1520    | 6" Dual Connection   | 32"-38"   | 92"L x 32.5"W x 32"-38"H |