



Covered Dissection Table with Lid Assist

HB Series

* Specifications of modified models may vary from the model shown.

Specifications

Standard Features

- Hinged covers with lid assist to safely and easily open and close
- All stainless steel construction
- Recessed, non-removable top
- Level cart surface with sloped recess
- Includes drain with PVC ball valve and drain tube
- Four heavy duty 8" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

Performance

- Hoods cover cadaver approximately up to 250 lbs.
- Weight capacity: Approximately 350 lbs

Optional Features

- Bridge style leg rests
- Bridge style book stand
- Removable bucket holder
- Stainless steel bucket
- Post style leg rests
- Post style book holder
- Removable recessed 48" shelf
- Lower fixed drawers



		Features		
	Model No.	Ventilation	Elevation	Dimensions
Product	HB1100			86"L x 32.5"W x 36"H
	HB1410	6" Single Connection		86"L x 32.5"W x 36"H
	HB1420	6" Dual Connection		86"L x 32.5"W x 36"H
	HB1510	6" Single Connection	32"-38"	92"L x 32.5"W x 32"-38"H
	HB1520	6" Dual Connection	32"-38"	92"L x 32.5"W x 32"-38"H