

Ensemble Dissection Board, Light Blue




AX066-048



AX066-060



AX066-072



AX066-096

SPECIFICATIONS

Construction

- Made to the depth of the grossing station working area
- The Ensemble Dissection Board is the width of one half of the grossing station work surface
- AX066-096 includes (2) dissection boards

Dimensions

- 18.5" D x 0.8" H
- AX066-048: 18" L
- AX066-060: 27" L
- AX066-072: 38" L
- AX066-096: 21" L

Capacity

- 40 lbs

Materials

- HDPE, High-Density polyethylene, light blue in color

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

Variations

Model No.	Description	Quantity	Dimensions
AX066-048	Dissection Board, Light Blue, 48" Units	1 Board	18" L x 18.5" D x 0.8" H
AX066-060	Dissection Board, Light Blue, 60" Units	1 Board	27" L x 18.5" D x 0.8" H
AX066-072	Dissection Board, Light Blue, 72" Units	1 Board	38" L x 18.5" D x 0.8" H
AX066-096	Dissection Board, Light Blue, 96" Units	2 Boards	21" L x 18.5" D x 0.8" H

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI