

Body Bag, USA Response, Envelope Zipper, Heavy Duty, PE, 4 Handles – BE172

Mopec has partnered with multiple vinyl and polyethylene suppliers to meet the demand for body bags. These new partnerships will not only enable us to meet the large quantities needed across the country, it will help keep US manufacturing companies operating and save American jobs. We are excited for these partnerships and the opportunity to take care of our valuable customers in the midst of this crisis.



SPECIFICATIONS

Standard Features

SIZE: Adult

DIMENSIONS: 36x96 inches

MATERIAL: 12 Mil Polyethylene

CONSTRUCTION: Heat sealed seams, leak proof.

WEIGHT CLASS: Heavy Duty

WEIGHT TESTED: 300lbs+

COLOR: Black

ZIPPER TYPE: Envelope (Curved) Zipper

CATEGORY: Disaster Bag, Human Remains Pouch

CHLORINE-FREE: Yes

HANDLES: No

THICKNESS: 12 Mil

SUITABLE FOR CREMATION: Yes

PRODUCTION: US



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071