

# Body Bag, USA Response, Envelope Zipper, Standard Duty, PVC, No Handles – BE173

Mopec has partnered with multiple vinyl and polyethylene suppliers to meet the demand for body bags. These new partnerships will not only enable us to meet the large quantities needed across the country, it will help keep US manufacturing companies operating and save American jobs. We are excited for these partnerships and the opportunity to take care of our valuable customers in the midst of this crisis.



## **SPECIFICATIONS**

### **Standard Features**

SIZE: Adult

DIMENSIONS: 36x93 inches

MATERIAL: PVC (Unsupported Vinyl)

CONSTRUCTION: Heat sealed seams, leak proof.

WEIGHT CLASS: Standard Duty

WEIGHT TESTED: 220lbs+

COLOR: Black/White/Varies

ZIPPER TYPE: Center (Straight) Zipper

CATEGORY: Disaster Bag, Human Remains Pouch

CHLORINE-FREE: Yes

HANDLES: No

THICKNESS: 7-10 Mil

SUITABLE FOR CREMATION: No

PRODUCTION: US



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071