

Mopec 810 Autopsy Saw

The Mopec 810 Autopsy saw remains as the industry standard in anatomic pathology. High speed oscillating action quickly cuts through bone with minimal damage to soft tissue. Used for removing the cranial cap, making linear cuts or sectioning small bone specimens.



SPECIFICATIONS

Construction

- Cast aluminum housing
- Arbor-mounted blades
- 10 foot cord with hospital grade plug
- Powerful 17,000 RPM motor produces 32,000 oscillations per minute
- Familiar lightweight & ergonomic design
- Slender handle for precise control

Materials

- Armature, field coil, cord, plugs, switches, bearings, carbon brushes, gaskets, cam, shaft, shaft link.

Dimensions/Weight

- Shipping Dimensions 6.5 x 18.5 x 5.5
- Weight with power cord 3.7 lbs
- Shipping weight: 4.2 lbs

Capacity

- 32,000 oscillations per minute

Utilities

- Voltage: 115V/60Hz/1PH

Included Options/Accessories

- Blades: BD101 Large Section Blade
- Additional Tools: Allen wrench

Mopec 810 Autopsy Saw Replacement Blades

Mopec Part #	Purpose/Description	Stryker Part #	Cut Edge x Cut Depth
BD101	Large Section Blade	1100	3.70" x 0.98"
BD104	Round Blade	1119	2.56" x 0.85"
BD105	Spinal Column Blade	1105	3.79" x 1.64"
BD108	Small Section Blade	1101	2.16" x 0.61"
BD110	Pediatric Blade	1106	0.79" x 1.62"
BD111	Deep Cutting Mastoid	1106A	0.80" x 1.62"



www.mopec.com
800.362.8491
21750 Coolidge Hwy.
Oak Park, MI 48237