


# Body Bag, Guardian Safe View, No Handles – BE175



## SPECIFICATIONS

### Standard Features

SIZE: Adult

DIMENSIONS: 36×96 inches (90×244 cm)

MATERIAL: Reinforced Coextruded Polyethylene

CONSTRUCTION: Heat sealed seams

WEIGHT CLASS: Standard Duty

WEIGHT TESTED: 300lbs

ZIPPER TYPE: Envelope (C-Shape) Zipper

CATEGORY: Body Bag

HANDLES: No Handles

THICKNESS: 10 Mil

ORIGIN: EU

COLOR: Black

### Design Notes

Unique optically clear interior panel allows viewing the remains contained in the bag with minimal exposure risk to personnel. Both panels are secured with envelope style heavy-duty zippers with dual zipper pulls for easy access.



[www.mopec.com](http://www.mopec.com)  
800.362.8491  
800 Tech Row  
Madison Heights, MI 48071