

Body Bag, US GOV/DOD Modified Design, Center Zipper, 18 Mil, 6 Handle – BE167-CZ

Using the design specifications for the Department of Defense (NATO) National Stock Number NSN-01-331-6244 Human Remains Pouch as their inspiration of design Mopec has secured exclusive production of 1,200 bags per day beginning Monday, April 13, 2020 and continuing through the end of the month.



SPECIFICATIONS

Standard Features

SIZE: Adult

DIMENSIONS: 36x96 inches (90x244 cm)

MATERIAL: PVC Supported Vinyl / Reinforced Vinyl Laminate

CONSTRUCTION: Heat sealed seams, leak proof.

WEIGHT CLASS: Heavy Duty

COLOR: Black

ZIPPER TYPE: Center (Straight) Zipper

ZIPPER DETAILS: Size #10

CATEGORY: Disaster Bag, Human Remains Pouch Heavy Duty

CHLORINE-FREE: No

HANDLES: 6 Handles w/ Bottom Support

THICKNESS: 17-18 Mil (minimum)

INNER LINER (BELOW BODY): No

CLEAR VIEWING PANEL: No

DESIGN: Designed to be similar to Dept of Defense (DOD) Specification NSN: 9930-01-331-6244

PRODUCTION: MEXICO & US

Deviations from Specification

Designed to be similar to Dept of Defense (DOD) Specification NSN: 9930-01-331-6244 in strength and durability, the Mopec BE167-CZ contains the following deviations in order to improve production time during the height of the COVID-19 crisis. Riveted handles also available.

ORIGINAL SPECIFICATION

Envelope Zipper
6 Riveted Handles

EXPEDITED PRODUCTION DEVIATION

Center Zipper
6 Heat Sealed Handles, with bottom support



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071