

# Level 4 Isolation Gown - BE304

With the ongoing COVID-19 pandemic, and an increase in demand for PPE, Mopec recognized a need for gowns that actually provide healthcare workers with an adequate level of protection and comfort. These level 4 isolation gowns offer the highest level of protection available and are used for high risk circumstances. The lightweight, breathable design allows for all-day comfort. They offer pathogen resistance, protection against non-airborne diseases, and are an effective barrier for fluids.



## SPECIFICATIONS

### Construction

- Level 4 Disposable
- SMS + PE composite
- Sterile for one time use
- Available in two sizes

### Materials

- 75% SMS + 25% PE

### Dimensions

- Size L: 48" L x 25.5" W
- Size XL: 50" L x 27" W
- Sleeves: 24" - 24.5"

### Color

- Blue

### Design

- Full back, wrap-around style
- Tie closures at neck and waist
- Long sleeves with elastic, white cuffs

### Precautions

- Use immediately after opening
- Limited to one-time use. Do not reuse.
- Discard immediately after use
- Product valid for three years



[www.mopec.com](http://www.mopec.com)  
800.362.8491  
800 Tech Row  
Madison Heights, MI 48071