


Mopec Biopsy Bags

Standard & Premium Grade Biopsy Bags Available. Ensure unrestricted fluid movement around tissue with the all new line of Mopec Biopsy Bags. Easy to open and available in three convenient sizes, you'll find our nylon mesh biopsy bags are the optimal solution for protecting specimens during tissue processing.



SPECIFICATIONS

Construction

- Nylon mesh, 0.2mm mesh diameter
- All seams are heat-sealed

Features

- Reduce reagent carryover and pressure artifacts often found with sponges
- Rapidly remove specimen with convenient peel-open bag
- Reduces the loss of small biopsy specimens

Options

- Three Standard Grade Sizes
- Two Premium Grade Sizes

ACTUAL SIZE REPRESENTATION



SKU	Description	Grade	Dimensions	Quantity
BI001	Biopsy Bags, Small	Standard	30mm X 50mm	1000/box
BI002	Biopsy Bags, Medium	Standard	45mm X 75mm	200/box
BI003	Biopsy Bags, Large	Standard	75mm X 95mm	500/box
BI004	Biopsy Bags, Small	Premium	30mm X 50mm	1000/box
BI005	Biopsy Bags, Medium	Premium	45mm X 75mm	200/box



www.mopec.com
 800.362.8491
 800 Tech Row
 Madison Heights, MI 48071