


# Morgue Refrigerator (10 Body) - KB1700

Drop-in self-contained refrigeration system with digital thermometer and easy to adjust temperature controller. Heavy-duty chrome plated hardware for long lasting use. Unit ships knocked down and will require assembly.



## SPECIFICATIONS

### Construction

- Front: Stainless steel with a no. 4 finish
- Top, sides, back, interior: Stucco embossed galvanized
- Interior flooring: Smooth aluminum with plywood reinforced understructure
- Insulation: 4" foamed in place urethane
- (10) Three-piece telescopic fabricated of extruded aluminum with nylon rollers
- Overlap door with magnetic gasket

### Materials

- Stainless steel and stucco embossed aluminum exterior
- Heavy-duty chrome-plated hardware for long lasting use
- Stainless steel removable trays

### Dimensions

- Overall: 190.5" W x 97" D x 80.5" H
- Doors: 27" W x 22" H

### Utilities / Power

- (2) 115V / 60Hz / 1 Ph
- (2) 208-230V / 50-60Hz / 1 Ph available (may affect dimensions & features)

### Controls

- Vapor proof light with switch

### Performance

- Drop-in self-contained refrigeration system with digital thermometer and easy to adjust temperature controller
- Operating Temperature: 38°F / 3.3°C

### Installation

- Unit ships knocked down and will require assembly

### Options & Accessories

- Custom sizes for larger cases
- Full length interior roller racks
- HD plastic storage boards



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071