


Morgue Refrigerator (3 Body), 24" Doors - KF1301

Drop-in self-contained refrigeration system with digital thermometer and easy to adjust temperature controller. Heavy-duty chrome plated hardware for long lasting use. Unit ships fully assembled - simply plug unit into an outlet and the unit is operational.



SPECIFICATIONS

Construction

- Exterior: Embossed galvanized steel construction
- Interior: Stucco embossed aluminum
- Interior flooring: Smooth aluminum with plywood reinforced understructure
- Insulation: 4" foamed in place urethane
- Provided with (3) interior roller racks
- Overlap door with magnetic gasket

Materials

- Stainless steel and stucco embossed aluminum exterior
- Heavy-duty chrome-plated hardware for long lasting use
- HD plastic storage boards

Dimensions

- Overall: 35" W x 88" D x 81" H
- Doors: 24" W x 52.75" H

Utilities / Power

- 115V / 60Hz / 1 Ph
- 208-230V / 50-60Hz / 1 Ph available (may affect dimensions & features)

Controls

- Vapor proof light with switch

Performance

- Drop-in self-contained refrigeration system with digital thermometer and easy to adjust temperature controller
- Operating Temperature: 38°F / 3.3°C

Installation

- Unit ships fully assembled, simply plug unit into an outlet and the unit is operational

Options & Accessories

- Custom sizes for larger cases
- Stainless steel body trays
- Casters: (4) 6" diameter casters, two are swivel type with wheel brakes



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071