

# Covered Dissection Table with Lid Assist

Designed for anatomy labs and teaching facilities. Features a recessed top which eliminates the need for a body tray. Built with a slope, the table drains toward one end and into an optional drain bucket neatly suspended under the frame. Two heavy-duty hinged covers are provided to securely conceal the cadaver in a discrete manner. During the dissection, the hinged covers open safely to protect from shin injury, and lock neatly and securely under the table. A PVC stop valve is located under the drain to stop the flow of fluids while emptying the bucket.



## SPECIFICATIONS

### Construction

- Heavy-duty hinged covers with lid assist to safely and easily open and close, and protect from shin injury
- Recessed, non-removable top
- Level cart surface with sloped recess
- Includes drain with PVC ball valve and drain tube

### Materials

- All stainless steel construction
- (4) Heavy duty 8" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

### Dimensions

- 86" L x 32.5" W or 92" L x 32.5" W
- Sloped top height: 35" - 36"

### Capacity

- Approximately 350 lbs

### Utilities / Power

- Exhaust requirements: 220 - 340 CFM
- Static pressure: .34 - .67" WG

### Performance

- Hoods cover cadavers up to approximately 350 lbs

### Options & Accessories

- Bridge-style or Post-style leg rests
- Bridge-style or Post-style book holder
- Removable bucket holder
- Stainless steel bucket
- Removable recessed 48" shelf
- Lower fixed drawers
- Misting system

### Features

Model No.	Ventilation	Elevation	Dimensions
HB1400			86"L x 32.5"W x 36"H
HB1410	6" Single Connection		86"L x 32.5"W x 36"H
HB1420	6" Dual Connection		86"L x 32.5"W x 36"H
HB1500		32"-38"	92"L x 32.5"W x 32"-38"H
HB1510	6" Single Connection	32"-38"	92"L x 32.5"W x 32"-38"H
HB1520	6" Dual Connection	32"-38"	92"L x 32.5"W x 32"-38"H



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071