

Self Propelled, Covered Cadaver Carrier - JA800

The Mopec JA800 Self Propelled, Covered Cadaver Carrier is designed to ensure your safety and efficiency while transporting discreetly. The motorized drive system virtually eliminates manual pushing and can be instrumental in reducing injuries during transferring and transporting. A false frame and cadaver cover are included to provide discreet transport from the hospital room to the mortuary.



SPECIFICATIONS

Construction

- Sturdy, lightweight, removable top frame works in conjunction with flip-down side rails to conceal cadaver
- Motorized drive system virtually eliminates manual pushing
- Motorized drive handles activate the drive system with a gentle squeeze and push of the handles for effortless transport
- Large wheel mobility can reduce start-up force up to 50% and steering effort up to 60%

Materials

- Top tray: 14 gauge stainless steel
- Cover: launderable cloth
- Top frame: aluminum

Dimensions

- Side Rails Up: 38" W x 86" L
- Side Rails Down: 31" W x 86" L
- Height (Lowest): 23.25" H
- Height (Highest): 34.5" H
- Height (Overall): 44.5"

Controls

- Brake controls located at the foot and head end for quick and easy accessibility
- Hydraulic controls located on either side of the stretcher for quick accessibility and efficiency

Performance

- Weight Capacity: 700 lbs
- Brakes on all casters lock both rotation and swivel for superior security

Utilities / Power

- 8' cord and retractable cord reel for charging zoom drive batteries
- 120V or 240V available

Options & Accessories

- Replacement cover available
- Integrated scale



www.mopec.com
 800.362.8491
 800 Tech Row
 Madison Heights, MI 48071