

Two Body Rotating Dissection Table - HM100

Specifically designed for use in anatomy labs where space is a concern with high student traffic. The two levels provide storage and access to the appropriate cadaver. The trays are smoothly rotated 180° with the anodized aluminum mechanism. Each tray has a small drain pan or drain plug.



SPECIFICATIONS

Construction

- Stainless steel construction
- (4) locking casters
- Includes (2) fixed body trays with a small drain pan or drain plug
- Aluminum mechanism for a smooth 180 degree rotation

Materials

- Stainless steel construction with heavy-duty 304 tubing
- Anodized aluminum

Performance

- Two tiers provide storage and access to two cadavers at once
- Smoothly rotates from tier position
- Rotates a full 180 degrees
- Weight capacity: 300 lbs per tray
- Cart should be loaded, unloaded, and rotated by two operators at all times to ensure a safe and smooth rotation
- The weight of the cadavers must be within approximately 20 lbs of each other

Dimensions

- Table: 95" L x 33" W x 39.7" H
- Trays: 78" L x 26" W



www.mopec.com
800.362.8491
800 Tech Row
Madison Heights, MI 48071