

Cadaver Lift, Fork Style, Powered - JD450 Series

The JD450 Series Fork Style Cadaver Lifts can lift up to 1000 lbs and can be raised as high as 64" or 76". These powered lifts alleviate the need of manual cranking to raise and lower their forks with a corded pendant with simple up/down elevation controls.



SPECIFICATIONS

Construction

- Heavy-duty painted steel
- (2) 5" Front phenolic wheels (static)
- (2) 4" Rear phenolic swivel casters

Dimensions

- JD450:
 - 25.5" Width
 - 44.5" Depth
 - 80" Height
- JD451:
 - 25.5" Width
 - 44.5" Depth
 - 93" Height

Elevation

- JD450:
 - Minimum Height: Floor
 - Maximum Lift Height: 64"
- JD451:
 - Minimum Height: Floor
 - Maximum Lift Height: 76"

Performance

- Safety screen
- Foot-operated floor lock
- Illuminated battery status gauge
- Load Capacity: 1000 lbs

Controls

- Corded pendant for elevation control

Options & Accessories

- JD450: 64" lift height
- JD451: 76" lift height
- JD460: Adjustable forks



www.mopec.com
800.362.8491
800 Tech Row
Madison Heights, MI 48071