



 Mopec

Formaway

FORMALIN NEUTRALIZER PRODUCTS

Formaway Formalin Neutralizing Absorbent Pads

Mopec's Formaway Formalin Neutralizing Absorbent Pad (F-NAP) was designed to reduce formaldehyde vapors while dissecting tissue, transferring tissue from containers of formalin, and transporting or storing tissue in containers of formalin, by absorbing and neutralizing formalin. Available in miniature, small, medium, and large. Sold by the pack or by the case.



| SKU | Description | Dimensions | Quantity |
|-------|------------------------------------|------------|----------|
| BE092 | Formalin Neutralizing Pads, Mini | 3" x 4" | 100/pack |
| BE093 | Formalin Neutralizing Pads, Small | 8" x 11" | 25/pack |
| BE094 | Formalin Neutralizing Pads, Medium | 11" x 15" | 25/pack |
| BE095 | Formalin Neutralizing Pads, Large | 15" x 20" | 25/pack |

Formaway Formalin Neutralizing Powder



Mopec's Formaway Formalin Neutralizing Powder safely and easily neutralizes spent formalin waste to a non-hazardous and non-toxic solution in only 20 minutes. Contents are emptied directly into 1 gallon of spent 10% Formalin or 4% Glutaldehyde, reducing risk exposure to lab personnel. This unique formula also eliminates drain-clogging polymers. Available in a case of 12 or 32 bottles.

| SKU | Description | Size | Quantity |
|----------|------------------------------|---------------|----------|
| BG550-12 | Formalin Neutralizing Powder | 12 oz. bottle | 12/case |
| BG550-32 | Formalin Neutralizing Powder | 12 oz. bottle | 32/case |
| BG551 | Waste Collection Container | 1 gallon | each |
| BG552 | Waste Collection Container | 2.5 gallons | each |
| BG553 | Funnel with Tissue Screen | - | each |
| BG554 | Formalin Test Strips | - | 100/pack |

Formaway Green Liquid Formalin Neutralizer

Mopec's Formaway Green Liquid Formalin Neutralizer quickly and effectively neutralizes spent formalin for sewer disposal. Simply add 32 ounces of Formaway Green to each gallon of 10% neutral buffered formalin waste and allow to stand for at least 4 hours before disposing. Available in 1 gallon, 2.5 gallon, and 5 gallon cubes.



| SKU | Description | Capacity |
|---------|-----------------------------------|-----------------|
| BG560-1 | Formalin Green Liquid Neutralizer | 1 gallon cube |
| BE560-2 | Formalin Green Liquid Neutralizer | 2.5 gallon cube |
| BE560-5 | Formalin Green Liquid Neutralizer | 5 gallon cube |



SCAN THE QR CODE OR VISIT MOPEC.COM/FORMAWAY FOR MORE INFORMATION