


APM Heavy Duty Body Bag, Envelope Zipper, 16 Mil PVC, 6 Handles | Optional Colors and Sizes - BE801E

Mopec is proud to partner with American Professional Manufacturing to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Custom-made, lockable slider for easy opening and closing
- Dielectrically sealed seams
- Envelope Zipper
- 6 handles secured by heavy-duty steel rivets
- Handles straps run completely under the pouch for added strength and support

Materials

- 16 Mil PVC Laminated Polyester
- HD nylon zipper, 100% sewn with nylon thread

Dimensions

- 94" L x 36" W

Capacity

- 300+ lbs

Performance

- Material will not tear, even if punctured
- 10 year shelf life
- Leak resistant bags that are impervious to blood and other fluids
- Autopsies and x-rays can be performed within the pouch

Options & Accessories

- Sold by the case of 6
- Available in black, red, blue, and yellow
- All color options are also available in the following sizes with 8 handles:
 - Adult XL, 8 Handles (94" L x 48" W)
 - Adult 2XL, 8 Handles (94" L x 60" W)
 - Adult 3XL, 8 Handles (94" L x 72" W)
 - Adult 4XL, 8 Handles (94" L x 84" W)
- Also available with a center zipper. Please contact a Mopec representative for more information.

Options

| Part No. | Dimensions | Capacity | Zipper Style | Color |
|-----------|---------------|----------|--------------|--------|
| BE801E-BK | 94" L x 36" W | 300+ lbs | Envelope | Black |
| BE801E-R | 94" L x 36" W | 300+ lbs | Envelope | Red |
| BE801E-BL | 94" L x 36" W | 300+ lbs | Envelope | Blue |
| BE801E-Y | 94" L x 36" W | 300+ lbs | Envelope | Yellow |

