


Standard Duty Body Bag, Envelope Zipper, 8 Mil PVC | Optional Colors and Sizes - BE802E

Mopec is proud to partner with American Professional Manufacturing to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Custom-made, lockable slider for easy opening and closing
- Dielectrically sealed seams
- Envelope Zipper

Materials

- 8 Mil PVC
- HD nylon zipper, 100% sewn with nylon thread

Dimensions

- 94" L x 36" W (BE802E)
- 94" L x 48" W (BE802E-XL)

Capacity

- 200+ lbs

Performance

- Leak resistant bags that are impervious to blood and other fluids
- 10 year shelf life

Options & Accessories

- Sold by the case of 12
- Available in white (BE802E) or yellow (BE802E-Y)
- Available in adult (BE802E) and bariatric (BE802E-XL) sizes.



Options

Part No.	Dimensions	Capacity	Zipper Style	Color
BE802E	94" L x 36" W	200+ lbs	Envelope	White
BE802E-Y	94" L x 36" W	200+ lbs	Envelope	Yellow
BE802E-XL	94" L x 48" W	200+ lbs	Envelope	White
BE802E-Y-XL	94" L x 48" W	200+ lbs	Envelope	Yellow

