

# Intermediate Duty Body Bag, Envelope Zipper, CF Polyethylene, Optional Handles

Mopec is proud to partner with American Professional Manufacturing to bring you the best body bags available in the market.



## SPECIFICATIONS

### Construction

- Custom-made, lockable slider for easy opening and closing
- Dielectrically sealed seams
- Envelope Zipper
- Optional Handles
- Handles straps run completely under the pouch for added strength and support

### Materials

- Blue polyethylene coated polyethylene fibers
- HD nylon zipper, 100% sewn with nylon thread
- Chlorine free

### Dimensions

- 94" L x 36" W

### Capacity

- 300 lbs

### Performance

- Material will not tear, even if punctured
- 10 year shelf life
- Leak resistant bags that are impervious to blood and other fluids
- Recommended strongly for cremation
- Autopsies can be performed in pouch

### Options & Accessories

- Sold by the case of 6 (BE811) or 10 (BE812)
- Available with (BE811) or without handles (BE812)

### Options

Part No.	Dimensions	Handles	Zipper Style	Color
BE811	94" L x 36" W	6 Handles	Envelope	Blue
BE812	94" L x 36" W	None	Envelope	Blue

