

Adult Medium Duty Disaster Pouch, Envelope Zipper, HDPE Tarp - BE411-224-EZ

Mopec is proud to partner with Salam International to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Ends are high frequency welded
- Heavy-duty #7 zipper
- 2 lockable pulls to secure pouch
- Complies with OSHA regulation for containment of body fluids, blood, body fats, chemicals, and fluids
- Envelope Zipper

Materials

- HDPE specially woven coated tarpaulin
- Chlorine free

Dimensions


- 94" L x 36" W

Performance

- 2 way zipper so pouch can open and close at any point
- No sewing or welding on sides to increase strength and eliminate leaks
- 100% leakproof
- Recommended for cremation

Options & Accessories

- Sold by the case of 20
- Available in blue, white and black




SALAM INTERNATIONAL, INC.

AVAILABLE FROM
