


Adult Medium Duty Disaster Pouch, Envelope Zipper, HDPE Tarp, 6 Handles - BE411-234-EZ

Mopec is proud to partner with Salam International to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Ends are high frequency welded
- Heavy-duty #7 zipper
- 2 lockable pulls to secure pouch
- Complies with OSHA regulation for containment of body fluids, blood, body fats, chemicals, and fluids
- Envelope Zipper
- Color: blue
- 6 Handles

Materials

- HDPE specially woven coated tarpaulin
- 2" wide webbing for easy grabbing
- Chlorine free

Dimensions

- 94" L x 36" W

Performance

- 2 way zipper so pouch can open and close at any point
- No sewing or welding on sides to increase strength and eliminate leaks
- Handles secured by heavy-duty rivets
- Handle straps run under the pouch for additional strength and support
- 100% leakproof
- Recommended for cremation

Options & Accessories

- Sold by the case of 15

