

Adult Heavy Duty Disaster Pouch, Envelope Zipper, Environmentally Friendly - BE411-384-EZ

Mopec is proud to partner with Salam International to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Ends are high frequency welded
- Heavy-duty #7 zipper
- 2 lockable pulls to secure pouch
- Complies with OSHA regulation for containment of body fluids, blood, body fats, chemicals, and fluids
- Envelope Zipper
- Color: white

Materials

- Ethylene copolymer mix with 3 layers of C4, C6, and C8
- Chlorine-free

Dimensions

- 94" L x 36" W

Performance

- 2 way zipper so pouch can open and close at any point
- No sewing or welding on sides to increase strength and eliminate leaks
- 100% leakproof
- Recommended for cremation

Options & Accessories

- Sold by the case of 20




SALAM INTERNATIONAL, INC.

AVAILABLE FROM

