

Adult Economy Disaster Pouch, Center or Envelope Zipper, 5-6 mil Vinyl

Mopec is proud to partner with Salam International to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Ends are high frequency welded
- Longer, heavy-duty #5 zipper
- 2 large lockable pulls to secure pouch
- Color: white

Materials

- 5-6 mil vinyl

Dimensions

- 90" L x 36" W

Performance

- No sewing or welding on sides to increase strength
- 2 way zipper with 2 sliders so bag can open and close in any location
- Wider envelope opening (BE411-124-EZ)
- 100% leakproof

Options & Accessories

- Sold by the case of 15
- Available with center zipper (BE411-122-CZ) or envelope zipper (BE411-124-EZ)



Options				
Part No.	Dimensions	Capacity	Zipper Style	Color
BE411-122-CZ	90" L x 36" W	200+ lbs	Center	White
BE411-124-EZ	90" L x 36" W	200+ lbs	Envelope	White

