

# Hydraulic Cadaver Carrier - JA600

The JA600s powerful pedal activated hydraulic system allows tray to be easily and quickly elevated to the height of the bed, helping to minimize back strain for the user. The JA600s form fitting, reinforced opaque cover over a sturdy lightweight aluminum frame ensures the utmost in discretion and consideration for all hospital personnel. Special rugged bumpers provide heavy-duty impact resistance. The JA600s Total Lock wheel locking system makes it easy to navigate and to securely lock the wheels in place for cadaver transfer. Lock one wheel in parallel alignment for easy steering through crowded halls and elevators. Wider body platform makes dealing with heavier bodies easier, yet the JA600 cart easily passes through doors and hallways. Contoured perimeter edges ensure retention and control of fluids.



## SPECIFICATIONS

### Construction

- Twin pedestal hydraulics for easy, foot-activated elevation
- (4) heavy-duty 8" casters with total lock
- Bumper guards
- Removable, aluminum false frame with heavy-duty, blue, opaque cover (JZ005)

### Materials

- All stainless steel fixed top

### Dimensions

- 81" L x 30" W

### Controls

- Brake lever on both sides and one end
- Elevating pedals on both sides

### Elevation

- Carrier only: 21.25" to 34"
- With false frame: 37.25" to 50"

### Capacity

- Weight capacity: 500 lbs.



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071